

HOLTER CONNECT™

Wireless Holter and Extended Holter Monitoring
Single-use Wearable Biosensors and
Cloud-based ECG Analysis and Reporting Service

Introducing the Latest Wireless ECG Patch Technology



Single-use Holter and Extended Holter Monitor

- Continuous 2-channel ECG, HR, Event recording
- Clinically-validated for hospital and home use
- Lowest cost of care on market



High Patient Compliance

- Minimal disruption to daily routine. Wearable during sleep, exercise and shower
- Easy application. No skin abrasion or shaving required
- Lightweight, only 28g



Clinically Reliable

- In a recent comparison study (N=141), arrhythmias were detected in a greater number of patients by the Wearable biosensor compared to the conventional Holter. *

Clinical Experience Highlights



100%
of participants found the biosensor to be "very comfortable" or "comfortable" to wear during the 7-day wear study**

92%
of participants had P-waveforms that were in the good to very good category and 8% were acceptable**


97%
of ECG data from patients was usable for analysis**

92%
of participants said "I often don't notice I'm wearing it"***

100%
of Biosensors data was retrieved with no packet loss**

0%
of Biosensors failed during the 7-day wear study**

Faster Reporting Faster Patient Diagnosis



- 1** **36 hours service level turn-around on Holter reports.**
Generated by certified cardiac technician service.
- 2** **Anytime access from any location.**
Access interactive reports, and full-disclosure ECG waveforms.
- 3** **Report commenting capability.**
Full authorization controls including E-signature & audit trail.

Patient-centric Anywhere Enrollment

IN CLINIC



Clinician applies Bio-sensor to Patient and enrolls and verifies ECG hook-up via handy App.

Patient wears Biosensor for prescribed time.

AT HOME



Biosensor assigned to Patient on Holter Connect Web Portal and mailed to Patient.

Patient self-applies Biosensor and wears for prescribed time.

Patient returns Biosensor. Data is uploaded to Holter Connect Secure Web Portal.

ABM processes data and provides Holter Report within 36 hours.

Clinical team can view, comment and e-sign Patient report.

Reimbursement and Billing

Holter

93224 Global
93225 Hook-Up
93226 Technical
93227 Professional

Extended Holter

93241 Global
93242 Hook-Up
93243 Technical
93244 Professional

All-Inclusive Rate

- Biosensor
- Report analysis & generation
- Platform access for unlimited number of users
- Data storage
- Training
- 24/7 service support

Two Billing Options

1. Split Bill
2. Fee for Service

To learn more, contact your local Advance Biometrix Representative at info@abmpatch.com or 1.877.859.4965

Holter Connect is a trademark of LifeSignals, Inc.
All product names, brands, logos, and trademarks are the property of their respective owners.

Advance Biometrix does not make any claims;
for informational purpose only.

Doc: ABM_1-2022

ADVANCE
BIOMETRIX

3304 NE 58th Street,
Vancouver, WA 98663, USA
advancebiometrix.com